

ALL INDIA INSTITUTE OF MEDICAL SCIENCES, BATHINDA Mandi Dabwali Road, Bathinda (Pb.) – 151001

PROCUREMENT CELL

AIIMS Bathinda/Proc.Cell/2025/071

Date: 07/04/2025

Due Date: 17/04/2025

Notice Inviting Quotations (NIQ)

The following items are required for AIIMS, Mandi Dabwali Road, Bathinda - 151001:

| Sr | Product with Description | Qty | Specifications |
|----|----------------------------|-----|----------------|
| 1 | Ballon Frontal Sinus Kit | 05 | Annexure – A |
| 2 | Ballon Maxillary Sinus Kit | 05 | Annexure – A |

Kindly arrange to send your quotation giving the lowest rates per unit along with terms and conditions in sealed cover addressed to the Executive Director, AIIMS Bathinda, Mandi Dabwali Road, Bathinda – 151001, Punjab so as to reach this office before the last date of receiving Quotations i.e., on 17/04/2025, Time 17:00

The word "Quotations for Ballon Frontal Sinus Kit & Ballon Maxillary for the dept. of ENT at AllMS Bathinda, our reference number and the date should boldly be mentioned on the cover.

- **1.** While submitting the quotation following should invariably be mentioned:
 - a) Details of Specification and its lowest rate F.O.R. AIIMS Bathinda.
 - b) GST Registration Number (attach copy)
 - c) PAN No. (attach copy)
 - d) Clearly mentioned GST rates.
 - e) Tax at a concessional rate as applicable to the Educational Institutions
 - f) Completion from the date of receipt of the confirmed order.
- 2. The rates will be valid for one year from the date of issuance of the Work Order.
- 3. Payment will be made after the supply and inspection of the delivered items.
- 4. Under no circumstances unsealed quotations will be entertained in the office.
- 5. Quotations received after the due date shall not be considered.

OIC, Procurement Cell

<u>UNDERTAKING</u>

(To be returned with quotation)

We hereby undertake the following:

- 1. We will not sell the product (s) to other institutions, bodies and also in the market at rates less than those quoted by us to AIIMS Bathinda.
- 2. The rate of GST/Excise Duty mentioned in the quotation is in accordance with the provisions of the rules and the same is payable to the concerned authorities in respect of
- 3. The goods/Stores/articles offered shall be of the best quality and workmanship and their supply will be strictly in accordance with the technical specifications and particulars as detailed in the quotation.
- The information furnished by us in the quotation is true and correct to the best of our knowledge and belief.
- 5. We have read and understood the rules, regulations, terms and conditions and agree to abide by them.

Authorised Signatory (Seal)

Quotation Proforma

(To be returned with Quotation)

Following proforma should be filled in and duly signed by the firm and sent along with the quotation.

| Fi | rm's Reference | | | | Date | | | |
|-----|-------------------------------|---|-----|---------------------------|------|--|--|----------------------|
| Fi | m Registration o. (if any) | | | | PAN | | | |
| GS | ST No. | | | Address: | | - Allendaria de la companya de la co | Comment of the Commen | |
| Ph | one No. | , | | | | | | |
| Em | ail Id | *************************************** | | | | | | |
| M/s | : | 7 ; | | NIQ No. | | | | |
| | | | | Last Date o Submission | | | | |
| Sr. | Product with | Description | Qty | Unit Price | GST | To Un Pri | Total Value | Total Value in words |
| 1 | | | | | | | | , |
| 2 | | | | ,- | | | | 1 |
| 3 | | | | | - | | | |
| , | Total | | | | | : 2 | | |

I/we engage to supply the material(s) as required to your office and comply the following:

- i. Terms and conditions of the NIQ.
- ii. That the offer is valid for 90 days.
- iii. That we have not been debarred by any Government body/PSU/Autonomous Body
- iv. That the rates quoted are not higher than the rates quoted for same item to any Government body/PSU/Autonomous Body.
- v. That the quotation submitted by us is properly sealed and prepared so as to prevent any subsequent alteration and replacement.

| Dlace | * | | Authorized Signatory |
|-----------------|---|--|----------------------|
| Place: Date: | | | (seal) |
| Dale. | | | |

Terms and Conditions

- 1. The quotation must be in the form furnished by Procuring Entity and should be free from corrections/erasures. In case there is any unavoidable correction it should be properly attested. If not, the quotation will not be considered. Quotation written in pencil will not be considered.
- 2. All India Institute of Medical Sciences (AIIMS), Bathinda reserves the right to accept the offer by individual items and reject any or all quotations without assigning any reason thereof and does not bind itself to accept lowest quotations.
- 3. Manufacturer's name and country of origin of materials offered must be clearly specified. Please quote whether your organisation is large scale industry or small scale industry. If you have NSIC/MSE/MSI/DGS&D Certificate, please attach it to the quotation. Mention your registration details.
- 4. Complete details and ISI specification if any must accompany the quotation. Make/brand of the item shall be stated wherever applicable. If you have got any counter offer as suitable to the material required by us, the same may be shown separately.
- 5. Samples must be submitted where specified along with the quotations. Samples must be carefully packed, sealed and labelled clearly with enquiry number, subject and sender's name for easy identification. Rejected samples will be returned at your cost if insisted.
- 6. All supplies are subject to inspection and approval before acceptance. Manufacturer/supplier warranty certificates and manufacturer/Government approved lab test certificate shall be furnished along with the supply, wherever applicable.
- 7. AIIMS Bathinda reserves the right to modify the quantity specified in this enquiry.
- 8. The prices quoted should be firm till the supplies are completed. Please quote the rates in words and figures. Rates quoted should be free delivery at destination including all charges otherwise the quotation is likely to be rejected. Prices quoted for free delivery at destination will be given preference. If there is no indication regarding the FOR, in the quotation, then it will be considered as FOR destinations. Price quoted should be net and valid for a minimum period of 90 days from the date of opening of the quotation.
- 9. Delivery period required for supplying the material should be invariably specified in the quotation.
- 10. In case your quotation is accepted and order is placed on you, the supply against the order should be made within the period stipulated in the order. AIIMS Bathinda reserves the right to recover any loss sustained due to delayed delivery by way of penalty. Failure to supply the material within the stipulated period shall entitle Procuring Entity for the imposition of penalty without assigning any reasons @ 0.5% (half percent) of the total value of the item covered in order as penalty per day subject to a maximum of 5% (five percent) unless extension is obtained in writing from the office on valid ground before expiry of delivery period.
- 11. If the deliveries are not maintained and due to that account Procuring Entity is forced to buy the material at your risk and cost from elsewhere, the loss or damage that may be sustained there by will be recovered from the defaulting supplier.
- 12. Dispute clause: Any dispute relating to the enquiry shall be subject to the jurisdiction of the court at Bathinda only.
- 13. Our normal payment terms are 100% within 30 (thirty) days on receipt and acceptance of material at our site in good condition.

BALLOON SINUS DILATION SYSTEM - BSDS

| SI.No | BALLOON SINUS DILATION SYSTEM - BSDS SPECIFICATION | SIZE |
|-------|--|--|
| 01112 | | |
| 1 | Sinus Balloon Catheter- It should have a non-complaint Balloon with puncture resistance. It should have a working length of 19cm. It should have a flexible shaft. It should have marker bands for endoscopic & direct visualization. It should have separate port for inflation / deflation and a separate port for introduction of the guide wire. Balloon should reach its stated size at 7 Atmosphere and should withstand 12 atmospheric pressure. Balloon Catheter should be approved by credible authorities for further safety of patients. Authorities like, US FDA, CE Approval & Indian DCGI. | 5mm x 17 mm length - Paediatric Balloons |
| | | |
| 2 | Sinus Balloon Catheter- It should have a non-complaint Balloon with puncture resistance. It should have a working length of 19cm. It should have a flexible shaft. It should have marker bands for endoscopic & direct visualization. It should have separate port for inflation / deflation and a separate port for introduction of the guide wire. Balloon should reach its stated size at 7 Atmosphere and should withstand 12 atmospheric pressure. Balloon Catheter should be approved by credible authorities for further safety of patients. Authorities like, US FDA, CE Approval & Indian DCGI. | 6mm x 17 mm length - Narrow Sinus Balloons for Adult |
| 3 | Sinus Balloon Catheter- It should have a non-complaint Balloon with puncture resistance. It should have a working length of 19cm. It should have a flexible shaft. It should have marker bands for endoscopic & direct visualization. It should have separate port for inflation / deflation and a separate port for introduction of the guide wire. Balloon should reach its stated size at 7 Atmosphere and should withstand 12 atmospheric pressure. Balloon Catheter should be approved by credible authorities for further safety of patients. Authorities like, US FDA, CE Approval & Indian DCGI. | 7mm x 17 mm length - Adul Balloons |
| 4 | Sinus Guide Catheter It should have a non-malleable stainless steel shaft. It should have clear distal portion for endoscopic visualization of balloon Sinuplasty instrumentation and a coloured atraumatic Tip for endoscopic visualization. It should consist of a suction barb and hub for clearing thin secretions. It should have a lumen to introduce balloon catheter. It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinuplasty kit. It should be compact in size to access tight nasal anatomy. Gright datheter should be approved by credible authorities for further safety of the should be approved by credible authorities for further safety of the should be should be approved by credible authorities for further safety of the should be should be should be approved by credible authorities for further safety of the should be sho | CUPTA |

Dept. of ENT & Head Neck

| and the contract of the contra | |
|--|---|
| Sinus Guide Catheter- It should have a non-malleable stainless steel shaft. It should have clear distal portion for endoscopic visualization of balloon sinuplasty instrumentation and a coloured atraumatic Tip for endoscopic visualization. It should consist of a suction barb and hub for clearing thin secretions. It should have a lumen to introduce balloon catheter. It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinuplasty kit. It should be compact in size to access tight nasal anatomy. Guide Catheter should be approved by credible authorities for further safety of patients. Authority like, US FDA, CE Approval & Indian DCGI. | 30' for Sphenoid sinus, |
| Sinus Guide Catheter- It should have a non-malleable stainless steel shaft. It should have clear distal portion for endoscopic visualization of balloon sinuplasty instrumentation and a coloured atraumatic Tip for endoscopic visualization. It should consist of a suction barb and hub for clearing thin secretions. It should have a lumen to introduce balloon catheter. It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinuplasty kit. It should be compact in size to access tight nasal anatomy. Guide Catheter should be approved by credible authorities for further safety of patients. Authority like, US FDA, CE Approval & Indian DCGI. | 70' for Frontal sinus |
| Sinus Guide Catheter- It should have a non-malleable stainless steel shaft. It should have clear distal portion for endoscopic visualization of balloon sinuplasty instrumentation and a coloured atraumatic Tip for endoscopic visualization. It should consist of a suction barb and hub for clearing thin secretions. It should have a lumen to introduce balloon catheter. It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinuplasty kit. It should be compact in size to access tight nasal anatomy. Guide Catheter should be approved by credible authorities for further safety of patients. Authority like, US FDA, CE Approval & Indian DCGI. | 90' for Maxillary sinus |
| It should have a non-malleable stainless steel shaft. It should have clear distal portion for endoscopic visualization of balloon sinuplasty instrumentation and a coloured atraumatic Tip for endoscopic visualization. It should consist of a suction barb and hub for clearing thin secretions. It should have a lumen to introduce balloon catheter. It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinus dilation kit. It should be compact in size to access tight nasal anatomy. Guide Catheter should be approved by credible authorities for further safety of patients. Authority like, US FDA, CE Approval & Indian DCGI. WKASPEEP GUPT | 110' for Maxillary sinus |
| | It should have a non-malleable stainloss steel shaft. It should have clear distal portion for endoscopic visualization of balloon sinuplasty instrumentation and a coloured atraumatic Tip for endoscopic visualization. It should bave a lumen to introduce balloon calheter. It should have a lumen to introduce balloon calheter. It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinuplasty kit. It should be compated in size to access light nasal anatomy. Guide Catheter should be approved by credible authorities for further safety of patients. Authority like, US FDA, CE Approval & Indian DCGI. Sinus Guide Catheter It should have a non-malleable stainless steel shaft. It should have clear distal portion for endoscopic visualization of balloon sinuplasty instrumentation and a coloured atraumatic Tip for endoscopic visualization. It should bave a lumen to introduce balloon catheter. It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinuplasty kit. It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinuplasty kit. It should be compated in size to access tight nasal anatomy. Guide Catheter should be approved by credible authorities for further safety of patients. Authority like, US FDA, CE Approval & Indian DCGI. Sinus Guide Catheter— It should have a non-malleable stainless steel shaft. It should be compated in size to access tight nasal anatomy. Guide Catheter should be approved by credible authorities for further safety of patients. Authority like, US FDA, CE Approval & Indian DCGI. Sinus Guide Catheter— It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinuplasty kit. It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinuplasty kit. It should be compatible with guide handle to enable single handed operation of the scope and Balloon sinuplasty kit. It should be compat |

| particular description of the second second | Sinus Illumination catheter- | |
|---|---|--|
| 9 | It should provide transcutaneous illumination for precise confirmation of access to paranasal sinuses. It should have an angle for ease of search in tight nasal anatomy. It should have a Floppy atraumatic tip to avoid damage to delicate nasal structure. It should have direct torque control with 1:1 ratio to enable precise placement of guide tip in the nasal anatomy. It should have a working length of more than 70 cm. It should be constructed of sturdy material to withstand multiple sinus access. It should have two optical channels to ensure bright transcutaneous illumination. It should have compatible accessories for all major light sources. Illumination wire should be approved by credible authorities for further safety of patients. Authority like, US FDA, CE Approval & Indian DCGI. | Illumination for Safe Navigation wire |
| | | |
| 10 | Sinus guide handle - It should be compatible with different sinus guide catheter. It should be malleable to suit different users. It should have a hand tight grip surface for better control and stability. It should have a lumen to introduce balloon sinus instrumentation. Guide handle should be approved by credible authorities for further safety of patients. Authority like, CE Approval & Indian DCGI. | Hand holding device |
| | | |
| 11 | Balloon inflation device — It should have an ergonomic piston handle for safe and controlled inflation/deflation of sinus balloon catheter. It should have clear grading of different pressure levels. It should have high pressure connector tubing to connect to the balloon catheter. It should be compatible with different sinus balloon catheters. It should have clear high pressure cylinder barrel to indicate the level of saline/ fluid inside the inflation device. It should have a piston locking system to prevent accidental inflation/deflation of the balloon catheter. It should be approved by credible authorities for further safety of patients. Authority like, CE Approval & Indian DCGI. | Balloon Inflation device |
| 12 | Light guide cable with adaptor — It should be compatible with all leading light sources. It should be compatible with sinus illumination system. It should have minimal heat transmission not exceeding body temperature at the tip of the sinus illumination system. It should be durable for multiple uses. | Light guide Cable |
| L was | ORDER PLACING CODES - KIT WISE. | |
| 13 | Balloon Maxillary Sinus Kits. Balloon Frontal Sinus Kits. Balloon Sphenoid Sinus Kit. Each Kit consist of 1 Balloon, 1 Guide Catheter, 1 Latch, 1 Illumination Wire, GUIT | - :A |
| Additi Dept. of | Each Kit consist of 1 Balloon, 1 Guide Catheter, 1 Latch 1 Illumination Wife, GUT Wikhil Rolessor & Head Sant Professor & Head Neck Surger ENT & Head Neck Surger ENT & Head Neck Surger ENT & Head Neck Surger Illis Bathinda Oept. of ENT & Head Neck Surger Associate Professor Associate Professor Associate Professor Associate Professor Associate Professor Associate Professor Allin's Bathinda | purgery |